

C2012 Haematocrit & Microtube centrifuge

*“looking for a centrifuge
that does both”*

Speed 500-12000Rpm

Max Rcf 13500G

Timer 0-99Minutes

Dims HWD 23x23x35Cms

Weight 12Kilos

Accessories for BR5401

8003 Haematocrit reader
8014 Sealing compound capillary



Rotors Available



BR5401

max speed 12000 Rpm 13,500G
24 x Capillary tubes
Tube size maximum H75 x W1.6mm
And 12 x 1.5-2.0ml tube
Tube size maximum H40 x W11mm
Cap size dependant



BR5424

Adaptors available
for BR5401 & BR5424

max speed 12000 Rpm 13,500G
24 x 1.5-2.2ml tube, fixed angle rotor
Tube size maximum H40 x W11mm
Cap size dependant
8004/6 - 0.2-0.4ml
8004/8 - 0.5ml



BR5436

max speed 12000 Rpm 13,500G
36 x 0.5ml tube, fixed angle rotor
Tube size maximum H40 x W8mm
Cap size dependant



BR5448

max speed 12000 Rpm 13,500G
48 x 0.2ml tube, fixed angle rotor
Tube size maximum H40 x W6mm
Cap size dependant