

C2015 Micro high speed centrifuge

"speedy little mover"

Speed 500-15000Rpm

Max Rcf 22000G

Timer 0-99Minutes

Dims HWD 23x23x35Cms

Weight 12Kilos



Rotors Available



BR5424

max speed 15000 Rpm 22,000G
24 x 1.5-2.2ml tube, fixed angle rotor
Tube size maximum H40 x W11mm
Cap size dependant
Adaptors available see below



BR5436

max speed 15000 Rpm 22,000G
36 x 0.5ml tube, fixed angle rotor
Tube size maximum H40 x W8mm
Cap size dependant



BR5448

max speed 15000 Rpm 22,000G
48 x 0.2ml tube, fixed angle rotor
Tube size maximum H40 x W6mm
Cap size dependant



8004/6

ADAPTORS for BR5424 only
24 adaptors for 0.2 -0.4ml tubes
Tube size maximum H40 x W6mm
Cap size dependant

8004/8

24 adaptors for 0.5ml tubes
Tube size maximum H40 x W8mm
Cap size dependant