

C2041 1 Litre Centrifuge

"from the family of workhorses"

Speed 500-4000Rpm

Max Rcf 2600G

Timer 0-99Minutes

Dims HWD 28x40x50Cms

Weight 22Kilos



Rotors Available



BR5624

max speed 4000 Rpm 2,250G
24 x 15ml tube, fixed angle rotor
Tube size maximum H120 x W17 mm
Cap size dependant
8009 - 5-10ml
8009/S - 2-4ml

Adaptors available



BR5608

max speed 4000 Rpm 2,250G
8 x 50ml tube, fixed angle rotor
Tube size maximum H115 x W30mm
Cap size dependant
6006 - 25ml
6007 - 2-4ml

Adaptors available



BR5510

max speed 4000 Rpm 2,600G
Swing out rotor 4 x 250ml max
Requires tube adaptors as below

B5250

Buckets required set of 4 (shown with Rotor)
Sealed lids available B5419



TUBE ADAPTORS for BR5510

Range from 0.5ml to 250ml max

See Page 14/15 for full details