

# C2043 3 Litre work horse

*"Big is best"*

Speed 500-4000Rpm

Max Rcf 3600G

Timer 0-99Minutes

Dims HWD 40x60x63Cms

Weight 43Kilos



## Rotors Available



### BR5553

max speed 4000 Rpm 3,600G  
4 x 750ml max swing out rotor  
Requires B3000 buckets as below



### B3000

max speed 4000 Rpm 3,600G  
Buckets for above set of 4  
Sealed lids available B5319

### ADAPTORS

for B3000

See Page 14/15 for full details



### B5200

max speed 4000 Rpm 2,800G  
Microtiter plate buckets  
Capacity 8 total stacked  
[requires BR5553 rotor - as shown]