

CR2000 Multi purpose & versatile rotor recognised centrifuge

"Top selling all rounder"

- Speed** 500-12000Rpm
- Max Rcf** 13500G
- Timer** 0-99Minutes
- Dims HWD** 27x33x46Cms
- Weight** 16.5Kilos



Rotors Available



BR5505

max speed 4000 Rpm 2,250G

4 x 100ml max swing out rotor with rings
Requires tube adaptors as below
Sealed lids available B5419



TUBE ADAPTORS for BR5505

Range from 0.5ml to 100ml max

See Page 14/15 for full details



BR5608

max speed 4000 Rpm 2,250G

8 x 50ml tube, fixed angle rotor
Tube size maximum H115 x W30mm
Cap size dependant

Adaptors available

6006 - 25ml
6007 - 15ml



BR5624

max speed 4000 Rpm 2,250G

24 x 15ml tube, fixed angle rotor
Tube size maximum H120 x W17 mm
Cap size dependant

Adaptors available

8009 - 5-10ml
8009/S - 2-4ml

Rotors - fixed angle



BR5401

max speed 12000 Rpm 13,500G
24 x Capillary tubes
Tube size maximum H75 x W1.6mm
Also 12 x 1.5-2.0ml tube
Tube size maximum H40 x W11mm
Cap size dependant.



BR5424

max speed 12000 Rpm 13,500G
24 x 1.5-2.2ml tube, fixed angle rotor
Tube size maximum H40 x W11mm
Cap size dependant.

Adaptors available
for BR5401 & BR5424

8004/6 - 0.2-0.4ml
8004/8 - 0.5ml



BR5436

max speed 12000 Rpm 13,500G
36 x 0.5ml tube, fixed angle rotor
Tube size maximum H40 x W8mm
Cap size dependant.



BR5448

max speed 12000 Rpm 13,500G
48 x 0.2ml tube, fixed angle rotor
Tube size maximum H40 x W6mm
Cap size dependant.



BR5412

max speed 6000 Rpm 3,800G
12 x 15ml tube, fixed angle rotor
Tube size maximum H120 x W17 mm
Cap size dependant.

Adaptors available

8009 - 5-10ml
8009/S - 2-4ml



BR5406

max speed 6000 Rpm 3,800G
6 x 50ml Tube, fixed angle rotor
Tube size maximum H115 x W30mm

Adaptors available

6006 - 25ml
6007 - 15ml