

## **K2015** *“Micro with the highest functionality & specifications”* **Ambient Model**

- Speed** 500-15000Rpm
- Max Rcf** 22000G
- Timer** 0-99Minutes + hold
- Dims HWD** 27x33x46Cms
- Weight** 16.5Kilos



## **K2015R** *“The micro chill...”* **Refrigerated Model**

- Speed** 500-15000Rpm
- Max Rcf** 22000G
- Timer** 0-99Minutes + hold
- Dims HWD** 35x40x62Cms
- Weight** 45Kilos
- Temp** -9 to +40C PID Controlled



## **Rotors - fixed angle** available to both



**BRK5424**

Adaptors available

max speed 15000 Rpm 22,000G  
24 x 1.5-2.2ml tube, fixed angle rotor  
Tube size maximum H40 x W11mm  
Cap size dependant  
8004/6 - 0.2-0.4ml  
8004/8 - 0.5ml



**BRK5436**

max speed 15000 Rpm 22,000G  
36 x 0.5ml tube, fixed angle rotor  
Tube size maximum H40 x W8mm  
Cap size dependant



**BRK5448**

max speed 15000 Rpm 22,000G  
48 x 0.2ml tube, fixed angle rotor  
Tube size maximum H40 x W6mm  
Cap size dependant



**BR5608**

Adaptors available

max speed 4000 Rpm 2,250G  
8 x 50ml tube, fixed angle rotor  
Tube size maximum H115 x W30mm  
Cap size dependant  
6006 - 25ml  
6007 - 15ml



**BR5624**

Adaptors available

max speed 4000 Rpm 2,250G  
24 x 15ml tube, fixed angle rotor  
Tube size maximum H120 x W17 mm  
Cap size dependant  
8009 - 5-10ml  
8009/S - 2-4ml