

# K2045 *"Designed with your every need in mind"*

Speed 500-12000Rpm

Max Rcf 13500G

Timer 0-99Minutes  
+ hold

Dims HWD 27x33x46Cms

Weight 17Kilos



## Rotors - swing out




### BR5505

max speed 4000 Rpm 2,250G  
4 x 100ml max swing out rotor  
with rings  
Requires tube adaptors as below  
Sealed lids available B5419



### TUBE ADAPTORS for BR5505

Range from 0.5ml to 100ml max  
See Page 14/15 for full details



### 8003

Reader for PCV values



### 8014

Capillary tube Sealant

Capillary tubes also available



### **BRK5401**

#### **Rotors - fixed angle**

max speed 12000 Rpm 13,500G  
24 x Capillary tubes  
Tube size maximum H75 x W1.6mm  
Also 12 x 1.5-2.0ml tube  
Tube size maximum H40 x W11mm  
Cap size dependant.



### **BRK5424**

max speed 12000 Rpm 13,500G  
24 x 1.5-2.2ml tube, fixed angle rotor  
Tube size maximum H40 x W11mm  
Cap size dependant.



### **BRK5436**

max speed 12000 Rpm 13,500G  
36 x 0.5ml tube, fixed angle rotor  
Tube size maximum H40 x W8mm  
Cap size dependant.



### **BRK5448**

max speed 12000 Rpm 13,500G  
48 x 0.2ml tube, fixed angle rotor  
Tube size maximum H40 x W6mm  
Cap size dependant.



### **BRK5624**

max speed 4000 Rpm 2,250G  
24 x 15ml tube, fixed angle rotor  
Tube size maximum H120 x W17 mm  
Cap size dependant  
8009 - 5-10ml  
8009/S - 2-4ml

Adaptors available



### **BRK5608**

max speed 4000 Rpm 2,250G  
8 x 50ml tube, fixed angle rotor  
Tube size maximum H115 x W30mm  
Cap size dependant  
6006 - 25ml  
6007 - 15ml

Adaptors available