

Prices remain until the stock is sold. No further discounts.
Normal 12 months warranty applies for all items.

NSW phone: (02) 9729 2988
VIC phone: (03) 9354 5788
QLD phone: 0434 678 246

Daihan WIG Series Laboratory Incubator Features

- Digital fuzzy control system for superior temperature accuracy
- Innovative easy-to-use controller built on door
- High quality LCD display with back light
- Simple control by jog shuttle switch
- Stainless steel chamber with corrosion resistant chrome shelves
- Over temperature protection, over current and circuit breaker



WIG-032 Laboratory Incubator - 32 Litres Capacity

Internal Dimensions	300x290x370 mm
External Dimensions	448x500x690 mm
Capacity	32 litres
Heater	150W

Stock: 2

Normally: \$1,450.00 plus GST

Now: \$1,300.00 plus GST

**Pricing does not include freight*



WIG-050 Laboratory Incubator - 50 Litres Capacity

Internal Dimensions	370x350x420 mm
External Dimensions	518x560x740 mm
Capacity	50 litres
Heater	200W

Stock: 2

Normally: \$1,600 plus GST

Now: \$1,440.00 plus GST

**Pricing does not include freight*



WIG-105 Laboratory Incubator - 105 Litres Capacity

Internal Dimensions	485x405x535 mm
External Dimensions	515x575x800 mm
Capacity	105 litres
Heater	200W

Stock: 2

Normally: \$1,879.00 plus GST

Now: \$1,695.00 plus GST

**Pricing does not include freight*



WIG-155 Laboratory Incubator - 155 Litres Capacity

Internal Dimensions	550x470x600 mm
External Dimensions	680x640x865 mm
Capacity	155 litres
Heater	300W

Stock: 2

Normally: \$2,150.00 plus GST

Now: \$1,935.00 plus GST

**Pricing does not include freight*