

Prices remain until the stock is sold. No further discounts.  
Normal 12 months warranty applies for all items.

Phone: (02) 9604 3911  
Email: [imports@thermoline.com.au](mailto:imports@thermoline.com.au)

## Daihan WIG Series Laboratory Incubator Features

- Digital fuzzy control system for superior temperature accuracy
- Adjustable temperature range from ambient to +70°C
- Innovative easy-to-use controller built on door
- High quality LCD display with back light
- Simple control by jog shuttle switch
- Stainless steel chamber with corrosion resistant chrome shelves
- Over temperature protection, over current and circuit breaker



### WIG-105, Laboratory Incubator (105 Litres Capacity)

Internal Dimensions	485x405x535 mm
External Dimensions	515x575x800 mm
Capacity	105 litres
Heater	200W

Stock: 1

Normally: \$2,066.00 plus GST

**Now: \$1,695.00 plus GST**

*\*Pricing does not include freight*

For more incubators, please view our complete product range at [www.thermoline.com.au/products.html](http://www.thermoline.com.au/products.html)