



-40°C Chest Freezer

For general lab applications



Technical Specs

Capacity:	121 Litres
Temperature Range:	-20°C to -45°C (manual defrost)
Temperature Control:	Thermostat with digital display
Dimensions External:	725mmL x 690mmD x 890mmH
Dimensions Internal:	520mmL x 345mmD x 630/390mmH
Basket:	1 Hanging (half width)
Lid:	Hinged with interior light, counter balanced hinges and key lock
External Finish:	White Enamel
Internal Finish:	Aluminium
Electrical:	230-240V / 50Hz / 1.4A
Optional Alarms:	Digital control with hi+lo temperature alarm (add \$500.00) Optional power failure alarm (add \$300.00) Optional BMS outlet to connect to customers building monitoring system (add \$150.00)

Special pricing applies while stocks last

\$2950.00 plus GST

(Only 2 units left)

Thermoline NSW (Head Office)
Email: sales@thermoline.com.au
9 Tarlington Place, Smithfield. NSW. 2164.
Phone: (02) 9604 3911
Fax: (02) 9725 1706

Thermoline VIC (Sales Office)
Email: vicsales@thermoline.com.au
11 Irene Street, Coburg North, VIC 3058. Australia
Phone: (03) 9354 5788
Fax: (03) 9354 5799

Thermoline NSW (Imports Office)
Email: imports@thermoline.com.au
Unit 12/55 Newton Rd, Wetherill Park. NSW. 2164.
Phone: (02) 9729 2988
Fax: (02) 9729 3699

CSK Climatek QLD & NT (Sales Office)
Email: sales@cskclimatek.com.au
P.O Box 3228, Darra, Brisbane QLD 4076
Phone: (07) 3376 5277

Thermoline: www.thermoline.com.au

LabDirect: www.labdirect.com.au