



# LABORATORY FREEZERS



# SPARK PROOF LABORATORY FREEZERS

(-20°C TO -25°C)

---

(Image: TF-365-SP-1-SD)



## Economy Laboratory Freezers

Static Defrost Freezers with Digital Temperature Control

- Fitted with digital controller
- High and low alarms standard
- Manual defrost (no defrost cycle)
- Key lockable doors
- Easy installation - plug in and set temperature
- Internal spark proof design
- BMS fitted standard

**Note: Only the interior is spark proof  
under normal operation.**





## Maximise your frozen storage space with Thermoline's spark proof freezer range.

Thermoline's spark proof freezers are fitted with a digital controller, high and low alarms, and are supplied standard with manual defrost (manual defrost ensures the temperature remains stable during operation). Shelves are in a fixed position and not adjustable.

The interior of this freezer is classified as 'spark free' as it does not have any wiring (such as a thermostat or light) which could ignite any vapours from hazardous materials which may be stored inside.

### Product Features

External connection for BMS (in-house building management system)

Lockable doors

Added benefit of digital control with indication of internal conditions

### Safety Features

High and Low Alarms



(Image: TF-140-SP-1-SD)

Model	TF-140-SP-1-SD	TF-365-SP-1-SD	TF-600-SP-1-SD
Capacity (litres)	140	365	600
Type	Upright		
Temperature Range	-20 °C to -25°C		
Storage Shelves	2	6	6
External Dimensions (WxDxH mm)	600x620x800	600x600x1850	780x695x1890
Internal Dimensions (WxDxH mm)	480x471x620	480x471x1620	657x581x1660
External Material	White Powder Coated Steel		
Internal Material	White Plastic		
Electrical	230V, 1A	230V, 2A	330V, 3.3A
Weight (kg)	44	84	90

**Head Office Phone** (02) 9604 3911

**Head Office Fax** (02) 9725 1706

**Email** [hello@thermoline.com.au](mailto:hello@thermoline.com.au)

**Web** [www.thermoline.com.au](http://www.thermoline.com.au)